

# ScopicLine™

## Disposable laparoscopic instruments A State of the art new generation of Hemostats made of plastic

Another innovative line of unique, single use, patented instruments from **Kencap Ltd.** *ScopicLine™* offers a series of clamps with various shapes and heads, designed to be used for: **Gynecology, Abdominal Surgery, General floor procedures, Emergency rooms, Doctor's office and Clinics.**

The *ScopicLine™* is used whenever a long and narrow instrument is needed for a stronger grip allowing an open clear view to see the procedures or to fit into small places.

*ScopicLine™* can come in different sizes and different heads. The heads and length can be designed to fit your customer needs. *ScopicLine™* can be sterilized by any type of sterilization, Gas (ETO), Autoclaving or Gamma Radiation. *ScopicLine™* is excellent for kit packers, reduce cost and help fight cross-infection.

### Advantages:

- State-of-the art design.
- Lock for holding in position.
- Large finger loops.
- No metal parts, Latex free.
- One dimension line for a open clear view.
- Disposable for one time use.
- Economically priced.
- Can be sterilized by all methods.
- Helps fight cross-infection.
- Made of medical approved material.

Overall length 300mm (12") with a 200mm (8") working arm. One unit weighs :25gr(0.055pound)



### ANATOMIC KELLY HEMOSTAT STRAIGHT

- Cat. No. 6343** *ScopicLine™* clamp straight not sterile bulk packed.
- Cat. No. 6340** *ScopicLine™* clamp straight without locking mechanism.
- Cat. No. 6345** *ScopicLine™* clamp straight sterile.



### ANATOMIC KELLY HEMOSTAT CURVED

- Cat. No. 6443** *ScopicLine™* clamp curved not sterile bulk packed.
- Cat. No. 6440** *ScopicLine™* clamp curved without locking mechanism.
- Cat. No. 6445** *ScopicLine™* clamp curved sterile.



### SPONGE HOLDER HEMOSTAT

- Cat. No. 6543** *ScopicLine™* sponge clamp not sterile bulk packed.
- Cat. No. 6540** *ScopicLine™* sponge clamp without locking mechanism.
- Cat. No. 6545** *ScopicLine™* sponge clamp sterile.




### PALMER TENACULUM CLASPING HEMOSTAT

- Cat. No. 6243** *ScopicLine™* Tenaculum clamp not sterile bulk packed.
- Cat. No. 6240** *ScopicLine™* Tenaculum clamp without locking mechanism.
- Cat. No. 6245** *ScopicLine™* Tenaculum clamp sterile.


**Packed:** 100 each in bulk or 50 sterile per box -36cm.X22cm.X14cm. (14"X8.5"X5.5").  
10 boxes per master case 46cm.X75cm.X39cm. (18"X29.5"X15.5")

**PATENTED**

## EXTENDED SCISSORS



Littauer



Sharp/Blunt



Sharp/Sharp



Blunt/Blunt

**CAT. No. 5425** SCISSORS LITTAUER.

**CAT. No. 5429** SCISSORS BLUNT/ BLUNT.

**CAT. No. 5427** SCISSORS SHARP/ BLUNT.

**CAT. No. 5423** SCISSORS SHARP/ SHARP.

All scissors are 220mm 8.5" long , 200 units per case.  
(CASE SIZE 36cm. X 22cm. X14cm.) 14"X8.5"X5.5" 10 CASES PER MASTER  
CARTON 46cm.X75cm.X39cm.

**All scissors can be packed individually in sterile pack. 100 packs per case**



## GAUGE MEASUREMENT (Hystometer)

**Cat. No.5222** not sterile

**Cat. No.5225** sterile

## FULL KIT CAT. No. 7820 Sterile Disposable

Contains: **ScopicLine™** Tenaculum palmer clamp, Scissors blunt/ blunt, **ScopicLine™** sponge clamp and measurement gauge.

## CAT. No. 7821 Sterile Disposable

Contains: **ScopicLine™** Tenaculum palmer clamp, Scissors blunt/ blunt, **ScopicLine™** clamp straight and measurement gauge.

## BASIC KIT CAT. No. 7822 Sterile Disposable –

Contains: **ScopicLine™** Tenaculum palmer clamp & Scissors blunt/ blunt.

**GYNIUS®PUNCH** is a unique disposable Biopsy Punch; A cervical biopsy is a test in which tissue samples are taken from the cervix and examined for disease or other problems. The Gynius® Punch is very stable & strong and has a special design for easy removal of biopsy sample after cut.

Disposable for One Time Use No cross contamination, Easy removal of biopsy sample after cut,  
Exact strong bite for clean tissue extraction, Made of medical approved material,  
Economically priced, Can be sterilized using any type of sterilization.



Easy removal of biopsy

**CAT No 5350 GYNIUS PUNCH®** Biopsy punch BULK NOT STERILE 100 bulk in a box

**CAT No 5355 GYNIUS PUNCH®** Biopsy punch STERILE Each one in a peel-pack pouch 50 per box

(CASE SIZE 36cm. X 22cm. X14cm.) 14"X8.5"X5.5" 10 CASES PER MASTER CARTON 46cm.X75cm.X39cm.

**KENCAP LTD.** P. O. Box 124, YAHUD 56101 ISRAEL, [kencap@kencapLtd.com](mailto:kencap@kencapLtd.com), [www.kencapLtd.com](http://www.kencapLtd.com)

**KENYON INDUSTRIES, INC** RINGWOOD, NJ. 07456 USA, [kenyon@warwick.net](mailto:kenyon@warwick.net)

[www.kencapLtd.com](http://www.kencapLtd.com)

PATENTED