



A WIDE RANGE OF DISPOSABLE DENTAL INSTRUMENTS AND KITS THESE DISPOSABLE DENTAL PRODUCTS COMBINE STATE OF THE ART PLASTIC AND METAL MATERIALS AND ENGINEERING, GIVE THE USER A STRONG GRIP OF HIGH QUALITY, COST EFFICIENT DISPOSABLE DENTAL INSTRUMENT.

DENTAL MIRROR



Double Sided Mirror

High quality and lightweight front surface mirror. One piece handle and mirror. Poly film protector with a pull-off tab. Available also with a double-sided mirror. Not made of glass, which increases the safety of the patient and eliminate surface reflection depth loss.

- Light weight
- One piece handle and mirror
- Poly film protector with a pull-off tab
- Available also with a double sided mirror
- Not made of glass, increases the safety of the patient
- No surface reflection depth loss

CAT. No. **6010** **ONE SIDED MIRRORS GREEN WITH TAB**
100 UNITS PER BOX (0.7kg per case)

CAT. No. **6012** **SAME AS CAT. NO. 6010 BUT WHITE**
72 UNITS PER BOX (0.5kg per case)

CAT. No. **6015** **SAME AS CAT. NO. 6010 BUT DOUBLE SIDED MIRROR**
100 UNITS PER CASE (0.9kg per case)

160mm (61/4") long, 10 BOXES PER MASTER CASE
CASE SIZE: 17cm x 12.5cm x 12.5cm (6.5"x5"x5") 10 CASES PER
MASTER CARTON, MASTER SIZE: 63cm x 26cm x 18cm (25"x10"x7")

COTTON & DRESSING FORCEPS

HIGH QUALITY STAINLESS STEEL AND PLASTIC

"MERIAM" type plier forceps. Tips inserted with formed, thin, tapered points of rigid plastic. Rigid plastic points will not scratch, score or damage other materials. No cool down or condensation. Perfect alignment of jaws with stabilizing pins. Plastic head would not cause irritation with patient, as metal does. Can be sterilized in any type of sterilization .

CAT. No. **3020**
150mm (6") long, 100 FORCEPCS PER BOX,
CASE SIZE: 23cm x 35cm x 15 (9"x13.7"x6"), 10 BOXES PER MASTER CASE,
MASTER SIZE: 73cm x 47cm x 39cm (28.7"x18.5"x15.3")

DENTAL EXPLORER

This dental explorer offers a new standard of explorer that has excellent performance and high quality stainless steel sharp edge. The handle is made of rigid plastic and designed in a diamond shape for improved function. Available also double side probe and Periodontal with colored edge.



Double side probe

- CAT. No. **6033** **PROBE DIAMOND # 23 EDGE,**
140mm (5½") long 100 UNITS PER BOX (0.4kg per case)
- CAT. No. **6033S** **SAME AS 6033 BUT STERILE**
100 UNITS PER BOX (0.4kg)
- CAT. No. **6037** **PROBE DIAMOND ANGLE SHAPE**
140mm (5½ ") 100 UNITS PER BOX (0.4kg per case)
- CAT. No. **6038** **PROBE PERIODONTAL WITH COLORED EDGE**
145mm (5¾ ") long, 100 UNITS PER BOX (0.4kg per case)
- CAT. No. **6041** **DOUBLE SIDE PROBE # 23 & ANGLE PROBE**
150mm (6") long, 100 UNITS PER BOX (0.6kg per case)
- CAT. No. **6043** **DOUBLE SIDE PROBE PERIODONTAL & #23 EDGE**
150mm (6") long, 100 UNITS PER BOX (0.6kg per case)

BOX SIZE: 17cm x 7cm x 12.5cm (6.5" x 2.8" x 5")
10 BOXES PER MASTER CARTON
MASTER SIZE: 49cm x 26cm x 19cm (19.2" x 10" x 7.5")


DENTAL



Probe diamond # 23 edge



Angle probe



Periodontal probe

KITS

Disposable cost efficient kits that contain two or three instruments; Mirror, Explorer & Forceps. Necessary for a routine dental examination of patients. Packed in a durable, transparent plastic peel pack tray, provides easy identification of contents.

- CAT. No. **7110** **NON- STERILE (MIRROR-PROBE & FORCEPS)**
50 UNITS PER CASE (1.9kg per case)
- CAT. No. **7112** **STERILE. (MIRROR-PROBE & FORCEPS)**
50 UNITS PER CASE (1.9kg per case)
- CAT. No. **7116** **STERILE. (MIRROR & PROBE ONLY)**
50 UNITS PER CASE (1.5kg per case)
- CAT. No. **7117** **STERILE. (MIRROR , FORCEPS AND**
DOUBLE SIDE PROBE PERIODONTAL & #23 EDGE)
- CAT. No. **7113** **STERILE. (MIRROR , FORCEPS AND**
DOUBLE SIDE PROBE ANGLE & #23 EDGE)

CASE SIZE: 43cm x 22cm x 14 (17" x 8.5" x 5.5")
10 CASES PER MASTER CARTON (500 kits)
MASTER SIZE: 47cm x 75cm x 46cm (18" x 29.5" x 15.5")
All kits are packed 50 trays in case

