



SCOPICLINE™

DISPOSABLE LAPAROSCOPIC INSTRUMENTS. A STATE OF THE ART NEW GENERATION OF HEMOSTATS MADE OF PLASTIC.

ScoplcLine™ offers a series of clamps with various shapes and heads, designed to be used for: Gynecology, General floor procedures, Emergency rooms, Doctor's office and Clinics.

The ScoplcLine™ is used whenever a long and narrow instrument is needed for a stronger grip allowing an open clear view to see the procedures or to fit into small places.

ScoplcLine™ can come in different sizes and different heads. The heads and length can be designed to fit your customer needs. ScoplcLine™ can be sterilized by any type of sterilization, Gas (ETO), Autoclaving or Gamma Radiation.

ScoplcLine™ is excellent for kit packers, reduces cost and helps fight cross-infection.



Unique, single use, patented instruments from Kencap Ltd.

The ScoplcLine™ is designed in a one dimensional line which allows an open clear view to see the procedures or to fit into the speculum. The ScoplcLine™ gives a stronger grip than conventional clamps while providing a more pleasant experience for the patient.



ScopIcLine™ can come in different sizes and different heads. The heads and length can be designed to fit your customer needs. ScopIcLine™ is excellent for kit packers, reduces cost and helps fight cross-infection.

Overall length 300mm (12") with a 200mm (8") working arm.

One unit weights: 25g (0.055 pound)

CASE SIZE: 23cm × 35cm × 15cm (9"×13.7"×6")

10 CASES PER MASTER CARTON (500 units)

MASTER SIZE: 39cm × 47cm × 73cm (15.3" × 18.5" × 28.7")

PALMER TENACULUM CLASPING HEMOSTAT

- Cat. No. **6243** **SCOPICLINE™ TENACULUM NOT STERILE BULK**
50 PER CASE (2.7 kg per case)
- Cat. No. **6240** **SCOPICLINE™ TENACULUM WITHOUT LOCKING**
50 PER CASE (2.7 kg per case)
- Cat. No. **6245** **SCOPICLINE™ TENACULUM STERILE**
50 PER CASE (2.5 kg per case)



ANATOMIC KELLY HEMOSTAT CURVED

- Cat. No. **6443** **SCOPICLINE™ CURVED NOT STERILE BULK**
50 PER CASE (2.7 kg per case)
- Cat. No. **6440** **SCOPICLINE™ CURVED WITHOUT LOCKING**
50 PER CASE (2.7 kg per case)
- Cat. No. **6445** **SCOPICLINE™ CURVED STERILE**
50 PER CASE (2.5 kg per case)



ANATOMIC KELLY HEMOSTAT STRAIGHT

- Cat. No. **6343** **SCOPICLINE™ CLAMP STRAIGHT NOT STERILE BULK**
50 PER CASE (2.7 kg per case)
- Cat. No. **6340** **SCOPICLINE™ CLAMP STRAIGHT WITHOUT LOCKING**
50 PER CASE (2.7 kg per case)
- Cat. No. **6345** **SCOPICLINE™ CLAMP STRAIGHT STERILE**
50 PER CASE (2.5 kg per case)



SPONGE HOLDER HEMOSTAT

- Cat. No. **6543** **SCOPICLINE™ SPONGE CLAMP NOT STERILE BULK**
50 PER CASE (2.7 kg per case)
- Cat. No. **6540** **SCOPICLINE™ SPONGE CLAMP WITHOUT LOCKING**
50 PER CASE (2.7 kg per case)
- Cat. No. **6545** **SCOPICLINE™ SPONGE CLAMP STERILE**
50 PER CASE (2.5 kg per case)





IUD INSTRUMENTS

In order to insert or remove IUDs the Doctor for this process needs instruments such as: Palmer claspings forceps, Scissors, Anatomic Kelly forceps and measuring gauge. You are selling to professionals, you don't need to explain what removal of IUD is, so don't need this part. We are proud to provide unique disposable instruments and kits for the insertion and/or removal of the IUD.

Our quality products are made of medical grade plastic for comfort for the patient and low chance of cross infection.

- Kits for full procedure
- Look more patient friendly

Cat. No. **5200** **PALMER CLASPING HEMOSTAT (TENACULUM)**

260mm (10.4") long
100 PER CASE (2.4kg per case)

Cat. No. **5205** **PALMER CLASPING IN STERILE PACK**

100 PER CASE (3.0 kg per case)



Cat. No. **5300** **ANATOMIC KELLY HEMOSTAT (SPONGE HOLDER)**

260mm (10.4") long
100 PER CASE (2.4kg per case)

Cat. No. **5305** **PALMER CLASPING IN STERILE PACK**

100 PER CASE (2.9 kg per case)



HYSTEROMETER GAUGE MEASUREMENT


232 mm (9") long , can bend and stay in the same position, plastic without memory.

Cat. No. **5222** **GAUGE MEASUREMENT**

300 PER CASE. (1.0 kg per case)

Cat. No. **5225** **GAUGE MEASUREMENT STERILE PACK**

200 PER CASE (1.7 kg per case)



CASE SIZE: 35cm x 22cm x 15cm (13.7" x 8.5" x 6")

10 CASES PER MASTER CARTON

MASTER SIZE: 47cm x 73cm x 39cm (18.5" x 28.5" x 15.5")

BIOPSY PUNCH

GYNIUS PUNCH® for single use disposable with stainless steel head to assures 100% Sterility clean, fast and accurate cutting and reduces the dangers of contamination.

A built in patented release extractor for taking out sample from biopsy.

A cervical biopsy is a test in which tissue samples are taken from the cervix and examined for disease or other problems. The biopsy identifies the abnormality when a positive pap smear shows minor cell changes or abnormalities.

A single use disposable GYNIUS PUNCH® with the new stainless steel metallic head and rigid plastic handles to assure 100% Sterility which reduces the dangers of contamination as well as infection and pain . The stainless steel head assure clean, fast and accurate cutting with minimum pain to the patient. Environmentally friendly.

For easy removal of biopsy sample, open handle (after cut) all the way until patented pin will enter the biopsy disk and sample will be extracted out to the glass plate or testing tube

ADVANTAGES

- Disposable for One Time Use
- Easy removal of biopsy sample after cut
- Exact bite for clean tissue extraction
- Firm hold of biopsy sample for confident procedures
- 3 guides to insure clean cutting and extra strength
- No cross contamination
- Easy and comfortable handles
- Spring grip for optimal handling
- Swivel lock to ensure sample retention
- Made of medical approved material
- Stainless steel rivet hold parts together
- Economically priced

CAT. No. 5350 GYNIUS PUNCH® BIOPSY PUNCH NOT STERILE

*228mm (11") long, Actual Bite Size: 5mm
50 BULK IN A BOX (2.2 kg per case)*

CAT. No. 5355 GYNIUS PUNCH® BIOPSY PUNCH STERILE

Each one in a peel-pack pouch 50 PER BOX (3 kg per case)

CASE SIZE: 35cm × 22cm × 15cm (13.7" × 8.5" × 6")

10 CASES PER MASTER CARTON (500 units)

MASTER SIZE: 39cm × 47cm × 73cm (15.3" × 18.5" × 28.7")



Release extractor for taking out sample



Stainless steel punch



IUD SCISSORS

Used for cutting IUD thread and various other medical procedures. Same as short scissors (page 4) but long handle (220mm, 8½" long)



Littauer



Sharp/Blunt



Sharp/Sharp



Blunt/Blunt

- CAT. No. **5425** **SUTURE REMOVAL**
200 PER CASE (3.6 kg per case)
- CAT. No. **5423** **SHARP / SHARP BLADES**
200 PER CASE (3.6 kg per case)
- CAT. No. **5427** **SHARP / BLUNT BLADES**
200 PER CASE (3.6 kg per case)
- CAT. No. **5429** **BLUNT / BLUNT BLADES**
100 PER CASE (2.0 kg per case)

CASE SIZE: 23cm x 35cm x 15 (9" x 13.7" x 6")
 10 CASES PER MASTER CARTON (2,000 units)
 MASTER SIZE: 39cm x 47cm x 73cm (15.3" x 18.5" x 28.7")

SKIN BIOPSY PUNCH

Skin biopsy is a biopsy technique in which a skin lesion is removed and sent to the pathologist to render a microscopic diagnosis. Round shaped knife ranging in size of 5mm is standard.

Instrument is 110mm (4½") long.

- CAT. No. **5365** **BIOPSY PUNCH 5MM**
100 PER CASE (0.5 kg per case)
- CAT. No. **5365S** **BIOPSY PUNCH 5MM STERILE**
100 PER CASE (0.6 kg per case)
- CAT. No. **5364** **BIOPSY PUNCH 4MM**
100 PER CASE (0.5 kg per case)
- CAT. No. **5364S** **BIOPSY PUNCH 4MM STERILE**
100 PER CASE (0.6 kg per case)
- CAT. No. **5363** **BIOPSY PUNCH 3MM**
100 PER CASE (0.5 kg per case)
- CAT. No. **5363S** **BIOPSY PUNCH 3MM STERILE**
100 PER CASE (0.6 kg per case)

CASE SIZE 17cm x 7cm x 12.5cm (6.5" x 2.8" x 5")
 10 BOXES PER MASTER CARTON
 MASTER SIZE: 49cm x 26cm x 19cm (9.2" x 10" x 7.5")